# New Farma Press Digital 8

Modello/Model/Modelo: FO3001

**AIR COMPRESSION THERAPY SYSTEM**

**NEW AGE ITALIA SRL**

Via De Brozzi, 3 – 48022 Lugo (RA)

Tel. 0545/32019

[www.newageitalia.it](http://www.newageitalia.it) – [info@newageitalia.it](mailto:info@newageitalia.it)

# New Farma Press Digital 8

## Manuale d'uso

Modello: FO3001

## PRESSOTERAPIA

### NEW AGE ITALIA SRL

Via De Brozzi, 3 – 48022 Lugo (RA)

Tel. 0545/32019

[www.newageitalia.it](http://www.newageitalia.it) – [info@newageitalia.it](mailto:info@newageitalia.it)

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**This manual has been drawn up to ensure user safety and exclude harmful uses.  
Please read this manual carefully before using our device.**

This device is for medical use

*This document is the property of New Age Italia srl. All rights reserved. It is forbidden to copy and reproduce by any means, including the total or partial photocopy of the content without written permission of New Age Italia srl.*

## INDEX

1. PRESENTATION .....	26
2. INDICATIONS.....	26
3. CONTRAINDICATIONS.....	26
4. SIDE EFFECTS .....	26
5. INFORMATION ON SAFETY .....	27
6. TECHNICAL SPECIFICATIONS.....	29
7. Product description .....	30
7.1 Standard endowment.....	30
7.2 Led control panel.....	31
7.2.1 Modes description.....	32
8. INSTRUCTIONS.....	34
8.1 How to connect the accessories.....	34
9. CLEANING AND MAINTENANCE .....	36
9.1 Cleaning .....	36
9.2 Maintenance.....	36
10. TROUBLESHOOTING .....	37
11. STORAGE.....	37
12. DISPOSAL .....	37
13. SYMBOLS.....	38
14. INFORMATION ON ELECTROMAGNETIC COMPATIBILITY (EMC).....	39

## 1. PRESENTATION

The pressotherapy system is characterized by the use of special automatic equipment using the properties of air compression on the human body and work together as a unit.

The air pump is connected to the accessories by a series of tubes. The direction of the compression massage goes from the extremity of the limb to the center of the body, inflating the air chambers in sequence and then deflating like a cycle. The pressure can be adjusted to avoid any discomfort to the patient. The accessories operate under the action of the sensor and microprocessor.

The pressotherapy system, in the medical field, thanks to the compression from the distal point to proximal, contributes to improve the circulation of blood and lymph, prevent deep vein thrombosis and relieve lymphatic edema.

## 2. INDICATIONS

The product is indicated for use by medical professionals and patients at home, who are under medical supervision. It helps to reduce swelling, prevent thrombosis in the limbs, treatment of lymphedemas related to tumors of the trunk or arm such as: primary lymphedema, edema as a result of trauma and sports injuries, post-immobilization edema, venous insufficiency, lymphedema. The system can be used at home or in hospital on patients aged 18 years or older.

## 3. CONTRAINDICATIONS

1. Acute inflammatory skin diseases
2. Arrhythmia
3. Erysipelas
4. Deep thrombophlebitis
5. Pulmonary edema
6. Acute venous thrombosis
7. Unstable hypertension
8. Patients with pacemakers
9. Presence of arterial insufficiency

## 4. SIDE EFFECTS


- Discomfort if the output pressure is too strong
- Redness or itching
- Sensations of constriction and pain during treatment with pressure above 120mmHg

## 5. INFORMATION ON SAFETY

*It is important to read all warnings and precautions included in this manual because they are intended to preserve the safety of the user, prevent injury and avoid situations that could cause damage to the device.*

- Do not use this device if you have a cardiac pacemaker, implanted defibrillator or other metal or electronic device implanted. Such use could cause electric shock, burns, electrical interference or death.
- Not suitable for use in presence of a flammable liquid
- The device should not be used by patients incapable of understanding and will.

<p><b>Repair by unauthorized personnel is not allowed. When the operation of the device is unusual, please do not intervene but contact assistance.</b></p> 	
<p><b>Stop using in case of anomalies and pain.</b></p>	<p><b>For the first use, the treatment time should not exceed 20 minutes.</b> For the first use, consult your doctor and do not use more than 20 minutes. Then gradually extend the treatment time according to the doctor's advice (long-term use may have side effects).</p>
<p><b>Store the device away from infants and young children.</b></p> <p>May cause damage such as electric shock.</p>	<p><b>Attention! Keep away from liquids and blunt objects.</b></p>
<p><b>The equipment should not be placed in places where it is difficult to disconnect the power supply.</b></p>	<p><b>Only use the 220-240V power connection and properly disconnect the power supply at the end of treatment.</b></p> <p>Improper power connection may cause damage to the device.</p>
<p><b>Handle with dry hands when connecting or disconnecting the power supply.</b></p> <p>Wet hands can cause electric shock</p> 	<p><b>Power cable can not withstand excessive stress, can not be placed under tables or chairs.</b></p> <p>If broken can cause electric shock.</p> 

<p><b>Do not bend the power supply cable.</b></p> 	<p><b>Keep accessories away from heat sources such as stoves, cigarettes etc.</b> They can cause damage to accessories.</p> 
<p><b>Do not use oil, benzene, alcohol, gasoline and chemicals to clean accessories.</b>          Otherwise, the accessories may be affected.          Clean the accessories with a damp cloth.</p>	<p><b>Do not open the zipper when the unit is working.</b>          This can cause damage to the accessory. Always turn off the device before removing the accessories.</p>
<p>Always use on clothes and underwear. Do not use on the naked body.</p> 	

## 6. TECHNICAL SPECIFICATIONS

DESCRIPTION	PARAMETERS
Number of chambers	8
Dimension of the device	240×200×120mm
Weight (with accessories)	2.22kg
Control system	LED touch
Treatment time	10~60 minutes
Treatment pressure	30~250mmHg
% Pressure error	<15%
Cycle time	85s~133s,±1s
Treatment modes	6 modes
Chamber setting	settable
Input voltage	AC 220-240V, 50-60Hz
Power consumption	65VA
Noise level	≤55dB
Protection against electric shock	Class II, Type BF Applied part
Mode of operation	Continuous operation
Degree of waterproof	IP21
Operating environment	+10°C - +40°C, 10% - 95%RH 86kPa-106kPa
Storage environment	-40°C - +70°C, 10% -100%RH 50kPa-106kPa
Dimension accessories	<b>Leg cuff:</b> Lenght 95cm Circumference: ankle 40cm, knee 60cm, hip 72cm
	<b>Extension for leg cuff:</b> added lenght 7cm
	<b>Waist cuff:</b> height 36 cm Circumference (90-115) cm
	<b>Arm cuff:</b> Lenght 85 cm
	<b>Total Body:</b> Total lenght: 125cm Leg lenght: 90cm Waist circumference: 100cm Thigh circumference: 64cm Calf circumference: 56cm

## 7. Product description

### 7.1 Standard endowment

Depending on the purchased version, the standard endowment is:

- A. NEW FARMA PRESS DIGITAL 8: device, accessories (2 leg cuffs, 1 waist cuff, 1 arm cuff), connexion cables for accessories and splitter cable for leg cuffs (the splitter allows the use of two leg cuffs at the same time)
- B. NEW FARMA PRESS DIGITAL 8 PLUS: device, accessories (2 leg cuffs, 1 waist cuff, 1 arm cuff), connexion cables for accessories, splitter cable for leg cuffs (the splitter allows the use of two leg cuffs at the same time) and TOTAL BODY with connecting cables and splitter cable
- C. NEW FARMA PRESS DIGITAL 8 BEAUTY: Device and TOTAL BODY (with connecting cables and splitter cable)

## 7.2 Led control panel

- **Time indicator (min)** : displaying the remaining treatment time
- **Time setting** : press to modify treatment time. Time range: 10~60 minuti with 10 min step. Default setting: 10 minutes. When time is over, machine stops automatically.
- **Pressure indicator** : indicating the selected treatment pressure through display and LED bead bar.
- **Pressure setting**: Pressure range 30-250mmHg. Default setting: 120mmHg. Press or to increase/decrease 10mmHg at a time.
- **Setting of sections**: default setting is all chambers active. Press to disable the corresponding chamber. Press again to enable.
- **On/off button**: Press to turn on. The device is then in stand-by; press again to turn off.
- **Start button** : when the device is in stand-by, the button is green-lighted. Press , and the button is blue-lighted. Press again to pause and the button will switch back to green light.
- **Mode setting**: 6 modes/programs are available (A) (B) (C) (D) (E) (F); default setting is mode (A). Press to change mode/program. Once the mode is selected, the corresponding indicators lighted.

## 7.2.1 Modes description

### **Standard mode** (peristaltic):

The chamber  $O_1$  inflates to the set pressure value within max 18 seconds. If the pressure of the chamber reaches the preset value within 18 seconds, it maintains the pressure for the rest of the period and then switches to the inflation of the next chamber. If the chamber pressure does not reach the preset value within the first 18 seconds, the inflation will be extended for an additional maximum of 12 seconds. Within this period of 12 seconds, if the pressure of the chamber reaches the preset value, it then switches to the inflation of the next chamber. If the chamber pressure cannot reach the preset value at the end of the additional 12 seconds, it automatically switches to the next chamber.

After 3 seconds, the next chamber (for example, chamber  $O_2$ ) inflates while the previous one deflates (for example chamber  $O_1$ ). The cycle of the 8 chambers repeats after the deflation of the last chamber 8

### **Sequential mode:**

The chamber  $O_1$  inflates to the set pressure value within a maximum of 18 seconds. If the chamber pressure reaches the preset value within 18 seconds, maintain the pressure and then move on to the inflation of the next chamber. If the chamber pressure does not reach the preset value within the first 18 seconds, the inflation will be extended for an additional maximum of 12 seconds. Within this period of 12 seconds, if the pressure of the chamber reaches the preset value, it then switches to the inflation of the next chamber. If the chamber pressure cannot reach the preset value at the end of the next 12 seconds, it automatically switches to the next chamber.

Once the inflation of chamber 8 has been completed, the chambers all deflate simultaneously. The next cycle starts again 13 seconds after the start of deflation.

### **Double mode:**

Chambers  $O_1 O_2$  inflate as a pair synchronously to the set pressure value within a maximum of 36 seconds. If the pressure reaches the preset value within 36 seconds, it switches to the inflation of the next torque. If the pair pressure does not reach the preset value within the first 36 seconds, the inflation will be extended for a maximum of 24 seconds. Within this period of 24 seconds, if the pressure of the pair of chambers  $O_1 O_2$  reaches the preset value, then switches to the next pair ( $O_3 O_4$ ). If the pressure of the chambers  $O_1 O_2$  cannot reach the preset value at the end of an additional 24 seconds, switch to the inflation of the next group pair ( $O_3 O_4$ ).

After 3 seconds, the next pair ( $O_3 O_4$ ) inflates, while the previous pair deflates. The next cycle begins 13 seconds after the last pair is inflated (chambers 7 e 8)


**D** Whole mode:

The 8 chambers inflate simultaneously at the set pressure within a maximum of 72 seconds and then deflate simultaneously. If the pressure does not reach the preset value within the first 72 seconds, the inflation is extended for a maximum of a further 28 seconds. At the set pressure level, the chambers deflate simultaneously. If the pressure does not reach the set level within 28 seconds, the chambers deflate automatically. The next cycle begins 13 seconds after the deflation of all chambers.

**E** A + B combined mode: standard (peristaltic) mode followed by the sequential one

**F** C + D combined mode: double mode followed by the whole one

## 8. INSTRUCTIONS

### 8.1 How to connect the accessories

Wear the accessories on light clothes.

#### A. Connect accessory cable pistons to accessory valves



**PLEASE NOTE:** the accessory cable consists of 8 cables with 8 terminations (pistons).

- Leg cuffs: connect the longest termination to the upper chamber (groin, corresponding to the piston and dark grey valve) and the others in succession;
- Arm cuff: connect the longest termination to the upper chamber (armpit, corresponding to the piston and valve of dark grey color) and the others in succession;
- Waist cuff: connect the longest termination to the upper chamber (chest, corresponding to the piston and dark grey valve) and the others in succession.
- Total body: connect the longest termination to the upper chamber (abdomen, corresponding to the piston and valve of dark grey color) and the others in succession.

*The accessory cables are identified by the length of the tubes and the distance between one piston and another. Depending on the accessory, the distance between a valve in one chamber and the next changes and consequently also the distance between the pistons in the respective cable.*

#### B. Insert the connector of the accessory cable to the air outlet of the device

**PLEASE NOTE:** connect by matching the parts marked with a small triangle

**PLEASE NOTE:** to use the two legs at the same time, connect the two cables corresponding to the splitter cable, and connect the latter to the air outlet of the device.

### C. Wear the accessories

**PLEASE NOTE:** to avoid elbow/knee pain, do not bend the limb when wearing leg and arm cuff. Before starting treatment, take a comfortable position, sitting or lying depending on the type of application

**D. Turn on the device,** set duration, pressure, number of chambers and treatment modes. Press  to start the treatment (the standby state is indicated with green light. The operating state is indicated with blue light).

## 8.2 Operation instructions

**IMPORTANT:** Before using the device, open the filter cover under the machine, remove the screw and put the lid back on. The screw is only used to protect the motor of the device from damage during transport.

1. Place on a flat surface
2. Connect the power supply to the power outlet (220~240V)
3. Connect accessory cables as described above
4. Connect the accessories cables to the air outlet of the device as above
5. Wear the accessories

**CAUTION:** make sure not to bend the cables or wrap them around the body

6. Press  to turn on the device
7. Set the pressure with  and .

*PLEASE NOTE: for the first use, it is recommended to select a pressure not more than 120 or 150 mmhg; then adjust the level depending on the sensation*

8. Set the time with 

Time frame is 10~60 minuti. The increase/decrease is 10 minutes. When the time ends, the device turns off automatically.

9. Select the mode with 

10. Disable chambers with button        .

*PLEASE NOTE: as default all chambers are active and corresponding indicators lighted. If you want to turn off a chamber, press the corresponding button ; press again to enable. When the button is not lit, it means that the chamber is disabled; if it flashes instead, the chamber is working.*

11. To pause during treatment, press . Press it again to resume treatment.

CAUTION: do not pull down the zipper of accessories during treatment, it can cause damage. Make sure to turn off the device before removing the accessories.

12. After treatment, disconnect the cable from the device

CAUTION: handle the cable that connects the accessories to the device carefully. Una volta connessi i cavi accessori ai gambali, fascia e bracciale, si consiglia di non staccarli per evitare danni. Once the accessory cables are connected to the leg, arm and waist cuff, it is advisable not to disconnect them to avoid damage.

## 9. CLEANING AND MAINTENANCE

The following indications are necessary to ensure the smooth operation of the product and to ensure its long-term safety.

### 9.1 Cleaning

Do not clean the device during treatment, make sure the device is turned off.

Clean the device with a soft, slightly damp cloth. In case of accumulation of heavy dirt, you can use a mild detergent.

Make sure that the water does not enter the machine. In this case, use the device only after it is completely dry.

Do not expose the device to moisture, do not keep it under running water, do not immerse it in water or other liquids.

Do not use any organic solvent or aggressive reagent for cleaning accessories. Do not wash the accessories with water. Clean the surface of the accessories thoroughly with a damp cloth.

### 9.2 Maintenance

1. In case of malfunctioning, contact the distributor or the point of sale.

We will not accept claims on repairs made by unauthorized personnel.

2. The user shall not repair the device or accessories. Contact the distributor or point of sale.

3. Opening/tampering by unauthorized third parties terminates all warranty claims

Each product has passed systematic validation at the manufacturing stage. Performance is stable and calibration and validation are not required. If your product does not achieve the expected performance and the basic function is changed during its normal use, please contact the distributor or point of sale.

**ASSISTANCE CENTER: New Age Italia s.r.l.**

Via de Brozzi 3 – 48022 Lugo (RA)

Tel. +39 0545 32019

www.newageitalia.it – info@newageitalia.it

## 10. TROUBLESHOOTING

Issue	Points to check
The device is not working	<ul style="list-style-type: none"><li>• Proper connection plug-switch</li><li>• Switching on of the device</li><li>• Proper power supply (220-240V )</li><li>• Automatic shutdown once treatment is over</li></ul>
Chamber is not inflating or there are leaks and unusual sounds	<ul style="list-style-type: none"><li>• If connectors are in good conditions</li><li>• If connector is well connected to the device and pistons are well connected to the valves of the accessories</li><li>• If connector cable is pressed or bent</li></ul>

## 11. STORAGE




Do not cut the connecting cables and accessories. Keep up away from sharp objects. After use, put the accessories in the original packaging. Do not expose the device to direct sunlight, protect from dirt and moisture. Store the device in a cool, well ventilated place. Never place heavy objects on the device.

## 12. DISPOSAL

At the end of the product life cycle, do not throw this product into normal household garbage, but take it to a collection point for recycling electrical equipment. Waste of electrical and electronic equipment can have potentially harmful effects on the environment. Incorrect disposal may cause accumulation of harmful toxins in air, water and soil and may be harmful to human health. You are obliged to dispose of the device correctly. Consult your municipality or dealer for disposal information.

## 13. SYMBOLS

SYMBOL	DESCRIPTION
	Information on the manufacturer Foshan Hongfeng Co.,Ltd, China
	EC Authorized Representative: MedNet EC-REP GmbH, Borkstrasse 10, 48163 Muenster, Germany
	Manufacturing date
	CE marking, Certificate issued by NB No.:0598
	Type of protection against electric shock: Class II equipment
	Degree of protection against electric shock: Type BF applied part
	Use indoors, not intended outdoor use
	Refer to the user manual
<b>IP21</b>	<b>Number 2:</b> protected against solid foreign objects of $\Phi 12.5\text{mm}$ and greater. <b>Number 1:</b> protected against falling water drops vertically
	<p>Numero di serie SN / serial nr. SN/ numero de serie SN</p> <p><input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/></p> <p>anno/year/ano      mese      numero di serie/ month      serial number/ mes      numero de serie</p>
	Electrical devices are recyclable materials and should not be disposed of with household waste after their useful life! Help us to protect the environment, save resources. Take the device to the appropriate collection points. Please contact the organization that is responsible for waste disposal in your area if you have any questions.
	Up to 4 layers
	FRAGILE, handle carefully
	Keep away from liquids and repair from rain

## 14. INFORMATION ON ELECTROMAGNETIC COMPATIBILITY (EMC)

<b>Guidance and manufacture's declaration – electromagnetic emission</b>			
FO3001 is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.			
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>	
RF emissions CISPR 11	Group 1	Its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emission CISPR 11	Class B	FO3001 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Class A		
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies		
<b>Guidance and manufacture's declaration – electromagnetic immunity</b>			
FO3001 is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	±1kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p>&lt;5% UT (&gt;95% dip in UT) for 0.5 cycle</p> <p>40% UT (60% dip in UT) for 5 cycles</p> <p>70% UT (30% dip in UT) for 25 cycles</p> <p>&lt;5% UT (&gt;95% dip in UT) for 5 sec</p>	<p>&lt;5% UT (&gt;95% dip in UT) for 0.5 cycle</p> <p>40% UT (60% dip in UT) for 5 cycles</p> <p>70% UT (30% dip in UT) for 25 cycles</p> <p>&lt;5% UT (&gt;95% dip in UT) for 5 sec</p>	<p>Mains power quality should be that of a typical commercial or hospital environment. If the user of the DEVICE requires continued operation during power mains interruptions, it is recommended that the DEVICE be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8</p>	<p>10V/m</p>	<p>10V/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>

NOTE UT is the a.c. mains voltage prior to application of the test level.

**Guidance and manufacture's declaration – electromagnetic immunity**

FO3001 is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	10V/m & table 9	10V/m & table 9	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> <p><math>d = 1.167\sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2.333\sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1) At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered.

If the measured field strength in the location in which the Device is used exceeds the applicable RF compliance level above, the Device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Device.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [Vi] V/m.

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment (Table 9)						
Test freq. (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	TETRA 400	Pulse modulation <sup>b)</sup> 18Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM <sup>c)</sup> ±5kHz deviation 1kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation <sup>b)</sup> 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217Hz	0.2	0.3	9
5500						
5785						

If necessary to achieve the immunity test level, the distance between the transmitting antenna and the FO3001 may be reduced to 1m. The 1m test distance is permitted by IEC 61000-4-3.

a. For some services, only the uplink frequencies are included.

b. The carrier shall be modulated using a 50 % duty cycle square wave signal.

c. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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**NewAge Italia Srl**